

Antibodies Product Datasheet

Product Name:	Mnk2 Antibody
Catalog Number:	C21220
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Mnk2.
Concentration:	1 mg/ml
Gene Name:	MKNK2
Protein Name:	MAP kinase-interacting serine/threonine-protein kinase 2
Other Name:	MKNK2; GPRK7; MNK2; MAP kinase-interacting serine/threonine-protein kinase 2; MAP kinase signal-integrating kinase 2; MAPK signal-integrating kinase 2; Mnk2
Human Swiss Prot No.:	Q9HBH9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mnk2 Polyclonal Antibody detects endogenous levels of Mnk2 protein.
Storage and Stability:	-20°C/1 year