

Antibodies Product Datasheet

Product Name:	ApoER2 Antibody
Catalog Number:	C21241
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LRP8. AA range:451-500
Concentration:	1 mg/ml
Gene Name:	LRP8
Protein Name:	Low-density lipoprotein receptor-related protein 8
Other Name:	LRP8; APOER2; Low-density lipoprotein receptor-related protein 8; LRP-8; Apolipoprotein E receptor 2
Human Swiss Prot No.:	Q14114
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ApoER2 Polyclonal Antibody detects endogenous levels of ApoER2 protein.
Storage and Stability:	-20°C/1 year