

Antibodies Product Datasheet

Product Name:	ERI1 Antibody
Catalog Number:	C21254
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ERI1. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	ERI1
Protein Name:	3'-5' exoribonuclease 1
Other Name:	ERI1; 3'EXO; THEX1; 3'-5' exoribonuclease 1; 3'-5' exonuclease ERI1; Eri-1 homolog; Histone mRNA 3'-end-specific exoribonuclease; Histone mRNA 3'-exonuclease 1; Protein 3'hExo; HEXO
Human Swiss Prot No.:	Q8IV48
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERI1 Polyclonal Antibody detects endogenous levels of ERI1 protein.
Storage and Stability:	-20°C/1 year