

Antibodies Product Datasheet

Product Name:	CD79a Antibody
Catalog Number:	C21274
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD79A. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	CD79A
Protein Name:	B-cell antigen receptor complex-associated protein alpha chain
Other Name:	CD79A; IGA; MB1; B-cell antigen receptor complex-associated protein alpha chain; Ig-alpha; MB-1 membrane glycoprotein; Membrane-bound immunoglobulin-associated protein; Surface IgM-associated protein; CD antigen CD79a
Human Swiss Prot No.:	P11912
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD79a Polyclonal Antibody detects endogenous levels of CD79a protein.
Storage and Stability:	-20°C/1 year