

Antibodies Product Datasheet

Product Name:	CD92 Antibody
Catalog Number:	C21279
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC44A1. AA range:391-440
Concentration:	1 mg/ml
Gene Name:	SLC44A1
Protein Name:	Choline transporter-like protein 1
Other Name:	SLC44A1; CD92; CDW92; CTL1; Choline transporter-like protein 1; CDw92; Solute carrier family 44 member 1; CD antigen CD92
Human Swiss Prot No.:	Q8WWI5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD92 Polyclonal Antibody detects endogenous levels of CD92 protein.
Storage and Stability:	-20°C/1 year