

Antibodies Product Datasheet

Product Name:	CD159a Antibody
Catalog Number:	C21289
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KLRC1. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	KLRC1
Protein Name:	NKG2-A/NKG2-B type II integral membrane protein
Other Name:	KLRC1; NKG2A; NKG2-A/NKG2-B type II integral membrane protein; CD159 antigen-like family member A; NK cell receptor A; NKG2-A/B-activating NK receptor; CD antigen CD159a
Human Swiss Prot No.:	P26715
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD159a Polyclonal Antibody detects endogenous levels of CD159a protein.
Storage and Stability:	-20°C/1 year