

Antibodies Product Datasheet

Product Name:	CD205 Antibody
Catalog Number:	C21294
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CD205.
Concentration:	1 mg/ml
Gene Name:	LY75
Protein Name:	Lymphocyte antigen 75
Other Name:	LY75; CD205; CLEC13B; Lymphocyte antigen 75; Ly-75; C-type lectin domain family 13 member B; DEC-205; gp200-MR6; CD antigen CD205
Human Swiss Prot No.:	O60449
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD205 Polyclonal Antibody detects endogenous levels of CD205 protein.
Storage and Stability:	-20°C/1 year