

Antibodies Product Datasheet

Product Name:	CD264 Antibody
Catalog Number:	C21301
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNFRSF10D. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	TNFRSF10D
Protein Name:	Tumor necrosis factor receptor superfamily member 10D
Other Name:	TNFRSF10D; DCR2; TRAILR4; TRUNDD; Tumor necrosis factor receptor superfamily member 10D; Decoy receptor 2; DcR2; TNF-related apoptosis-inducing ligand receptor 4; TRAIL receptor 4; TRAIL-R4; TRAIL receptor with a truncated death domain; CD
Human Swiss Prot No.:	Q9UBN6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD264 Polyclonal Antibody detects endogenous levels of CD264 protein.
Storage and Stability:	-20°C/1 year