

Antibodies Product Datasheet

Product Name:	CD272 Antibody
Catalog Number:	C21303
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CD272.
Concentration:	1 mg/ml
Gene Name:	BTLA
Protein Name:	B- and T-lymphocyte attenuator
Other Name:	BTLA; B- and T-lymphocyte attenuator; B- and T-lymphocyte-associated protein; CD antigen CD272
Human Swiss Prot No.:	Q7Z6A9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD272 Polyclonal Antibody detects endogenous levels of CD272 protein.
Storage and Stability:	-20°C/1 year