

Antibodies Product Datasheet

Product Name:	Atg9a Antibody
Catalog Number:	C21323
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of Atg9a
Concentration:	1 mg/ml
Gene Name:	ATG9A APG9L1
Protein Name:	Atg9a
Other Name:	ATG9 autophagy related 9 homolog A (S. cerevisiae)
Human Swiss Prot No.:	Q7Z3C6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Atg9a Polyclonal Antibody detects endogenous levels of Atg9a
Storage and Stability:	-20°C/1 year