

Antibodies Product Datasheet

Product Name:	Atg16L2 Antibody
Catalog Number:	C21326
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATG16L2. AA range:241-290
Concentration:	1 mg/ml
Gene Name:	ATG16L2
Protein Name:	Autophagy-related protein 16-2
Other Name:	ATG16L2; WDR80; Autophagy-related protein 16-2; APG16-like 2; WD repeat-containing protein 80
Human Swiss Prot No.:	Q8NAA4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Atg16L2 Polyclonal Antibody detects endogenous levels of Atg16L2 protein.
Storage and Stability:	-20°C/1 year