

Antibodies Product Datasheet

Product Name:	GABARAP Antibody
Catalog Number:	C21329
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of GABARAP
Concentration:	1 mg/ml
Gene Name:	GABARAPL2 FLC3A GEF2
Protein Name:	GABARAP
Other Name:	GABA(A) receptor-associated protein-like 2
Human Swiss Prot No.:	P60520
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABARAP Polyclonal Antibody detects endogenous levels of GABARAP
Storage and Stability:	-20°C/1 year