

Antibodies Product Datasheet

Product Name:	IL-1R9 Antibody
Catalog Number:	C21349
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL1RAPL2. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	IL1RAPL2
Protein Name:	X-linked interleukin-1 receptor accessory protein-like 2
Other Name:	IL1RAPL2; IL1R9; X-linked interleukin-1 receptor accessory protein-like 2; IL-1 receptor accessory protein-like 2; IL-1-RAPL-2; IL-1RAPL-2; IL1RAPL-2; IL1RAPL-2-related protein; Interleukin-1 receptor 9; IL-1R-9; IL-1R9; Three immunoglobuli
Human Swiss Prot No.:	Q9NP60
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-1R9 Polyclonal Antibody detects endogenous levels of IL-1R9 protein.
Storage and Stability:	-20°C/1 year