

Antibodies Product Datasheet

Product Name:	ECM1 Antibody
Catalog Number:	C21374
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human ECM1.
Concentration:	1 mg/ml
Gene Name:	ECM1
Protein Name:	Extracellular matrix protein 1
Other Name:	ECM1; Extracellular matrix protein 1; Secretory component p85
Human Swiss Prot No.:	Q16610
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ECM1 Polyclonal Antibody detects endogenous levels of ECM1 protein.
Storage and Stability:	-20°C/1 year