

Antibodies Product Datasheet

Product Name:	Dsg1 Antibody
Catalog Number:	C21388
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DSG1. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	DSG1
Protein Name:	Desmoglein-1
Other Name:	DSG1; CDHF4; Desmoglein-1; Cadherin family member 4; Desmosomal glycoprotein 1; DG1; DGI; Pemphigus foliaceus antigen
Human Swiss Prot No.:	Q02413
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dsg1 Polyclonal Antibody detects endogenous levels of Dsg1 protein.
Storage and Stability:	-20°C/1 year