

Antibodies Product Datasheet

Product Name:	LHX1 Antibody
Catalog Number:	C21410
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human LHX1.
Concentration:	1 mg/ml
Gene Name:	LHX1
Protein Name:	LIM/homeobox protein Lhx1
Other Name:	LHX1; LIM-1; LIM1; LIM/homeobox protein Lhx1; LIM homeobox protein 1; Homeobox protein Lim-1; hLim-1
Human Swiss Prot No.:	P48742
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LHX1 Polyclonal Antibody detects endogenous levels of LHX1 protein.
Storage and Stability:	-20°C/1 year