

## Antibodies Product Datasheet

|                               |                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Contactin 1 Antibody                                                                                                  |
| <b>Catalog Number:</b>        | C21420                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                     |
| <b>Applications:</b>          | WB, ELISA                                                                                                             |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.                                   |
| <b>Modification:</b>          | Unmodified/Total                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from the N-terminal region of human Contactin 1.                                          |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                               |
| <b>Gene Name:</b>             | CNTN1                                                                                                                 |
| <b>Protein Name:</b>          | Contactin-1                                                                                                           |
| <b>Other Name:</b>            | CNTN1; Contactin-1; Glycoprotein gp135; Neural cell surface protein F3                                                |
| <b>Human Swiss Prot No.:</b>  | Q12860                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Specificity:</b>           | Contactin 1 Polyclonal Antibody detects endogenous levels of Contactin 1 protein.                                     |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                          |