

Antibodies Product Datasheet

Product Name:	EBP1 Antibody
Catalog Number:	C21421
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PA2G4. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	PA2G4
Protein Name:	Proliferation-associated protein 2G4
Other Name:	PA2G4; EBP1; Proliferation-associated protein 2G4; Cell cycle protein p38-2G4 homolog; hG4-1; ErbB3-binding protein 1
Human Swiss Prot No.:	Q9UQ80
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EBP1 Polyclonal Antibody detects endogenous levels of EBP1 protein.
Storage and Stability:	-20°C/1 year