

Antibodies Product Datasheet

Product Name:	Chymase Antibody
Catalog Number:	C21430
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Chymase.
Concentration:	1 mg/ml
Gene Name:	CMA1
Protein Name:	Chymase
Other Name:	CMA1; CYH; CYM; Chymase; Alpha-chymase; Mast cell protease I
Human Swiss Prot No.:	P23946
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Chymase Polyclonal Antibody detects endogenous levels of Chymase protein.
Storage and Stability:	-20°C/1 year