

Antibodies Product Datasheet

Product Name:	DBP Antibody
Catalog Number:	C21435
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from DBP . at AA range: 90-170
Concentration:	1 mg/ml
Gene Name:	GC
Protein Name:	Vitamin D-binding protein
Other Name:	GC; Vitamin D-binding protein; DBP; VDB; Gc-globulin; Group-specific component
Human Swiss Prot No.:	P02774
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DBP Polyclonal Antibody detects endogenous levels of DBP protein.
Storage and Stability:	-20°C/1 year