

Antibodies Product Datasheet

Product Name:	Fyb Antibody
Catalog Number:	C21444
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Fyb . at AA range: 620-700
Concentration:	1 mg/ml
Gene Name:	FYB
Protein Name:	FYN-binding protein
Other Name:	FYB; SLAP130; FYN-binding protein; Adhesion and degranulation promoting adaptor protein; ADAP; FYB-120/130; p120/p130; FYN-T-binding protein; SLAP-130; SLP-76-associated phosphoprotein
Human Swiss Prot No.:	O15117
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fyb Polyclonal Antibody detects endogenous levels of Fyb protein.
Storage and Stability:	-20°C/1 year