

Antibodies Product Datasheet

Product Name:	Glypican-2 Antibody
Catalog Number:	C21487
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Glypican-2
Concentration:	1 mg/ml
Gene Name:	GPC2
Protein Name:	Glypican-2
Other Name:	Glypican-2 [Cleaved into: Secreted glypican-2]
Human Swiss Prot No.:	Q8N158
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Glypican-2
Storage and Stability:	-20°C/1 year