

Antibodies Product Datasheet

Product Name:	Internexin- α Antibody
Catalog Number:	C21539
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Internexin- α .
Concentration:	1 mg/ml
Gene Name:	INA
Protein Name:	Alpha-internexin
Other Name:	INA; NEF5; Alpha-internexin; Alpha-Inx; 66 kDa neurofilament protein; NF-66; Neurofilament-66; Neurofilament 5
Human Swiss Prot No.:	Q16352
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Internexin- α Polyclonal Antibody detects endogenous levels of Internexin- α protein.
Storage and Stability:	-20°C/1 year