

Antibodies Product Datasheet

Product Name:	Dectin-2 Antibody
Catalog Number:	C21598
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Dectin-2.
Concentration:	1 mg/ml
Gene Name:	CLEC6A
Protein Name:	C-type lectin domain family 6 member A
Other Name:	CLEC6A; CLECSF10; DECTIN2; C-type lectin domain family 6 member A; C-type lectin superfamily member 10; Dendritic cell-associated C-type lectin 2; DC-associated C-type lectin 2; Dectin-2
Human Swiss Prot No.:	Q6EIG7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dectin-2 Polyclonal Antibody detects endogenous levels of Dectin-2 protein.
Storage and Stability:	-20°C/1 year