

Antibodies Product Datasheet

Product Name:	HPLN1 Antibody
Catalog Number:	C21609
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HPLN1
Concentration:	1 mg/ml
Gene Name:	HAPLN1 CRTL1
Protein Name:	HPLN1
Other Name:	Hyaluronan and proteoglycan link protein 1 (Cartilage-linking protein 1) (Cartilage-link protein) (Proteoglycan link protein)
Human Swiss Prot No.:	P10915
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HPLN1
Storage and Stability:	-20°C/1 year