

Antibodies Product Datasheet

Product Name:	CD156c Antibody
Catalog Number:	C21625
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CD156c.
Concentration:	1 mg/ml
Gene Name:	ADAM10
Protein Name:	Disintegrin and metalloproteinase domain-containing protein 10
Other Name:	ADAM10; KUZ; MADM; Disintegrin and metalloproteinase domain-containing protein 10; ADAM 10; CDw156; Kuzbanian protein homolog; Mammalian disintegrin-metalloprotease; CD156c
Human Swiss Prot No.:	O14672
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD156c Polyclonal Antibody detects endogenous levels of CD156c protein.
Storage and Stability:	-20°C/1 year