

Antibodies Product Datasheet

Product Name:	FDC-SP Antibody
Catalog Number:	C21627
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human FDC-SP.
Concentration:	1 mg/ml
Gene Name:	FDCSP
Protein Name:	Follicular dendritic cell secreted peptide
Other Name:	FDCSP; C4orf7; Follicular dendritic cell secreted peptide; FDC secreted protein; FDC-SP
Human Swiss Prot No.:	Q8NFU4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FDC-SP Polyclonal Antibody detects endogenous levels of FDC-SP protein.
Storage and Stability:	-20°C/1 year