

Antibodies Product Datasheet

Product Name:	Hepassocin Antibody
Catalog Number:	C21633
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Hepassocin.
Concentration:	1 mg/ml
Gene Name:	FGL1
Protein Name:	Fibrinogen-like protein 1
Other Name:	FGL1; HFREP1; Fibrinogen-like protein 1; HP-041; Hepassocin; Hepatocyte-derived fibrinogen-related protein 1; HFREP-1; Liver fibrinogen-related protein 1; LFIRE-1
Human Swiss Prot No.:	Q08830
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Hepassocin Polyclonal Antibody detects endogenous levels of Hepassocin protein.
Storage and Stability:	-20°C/1 year