

Antibodies Product Datasheet

Product Name:	MAD2 Antibody
Catalog Number:	C21657
Reactivity:	Human, Rat, Mouse
Applications:	WB, IP, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IP 1:200-500 ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAD2L1. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	MAD2L1
Protein Name:	Mitotic spindle assembly checkpoint protein MAD2A
Other Name:	MAD2L1; MAD2; Mitotic spindle assembly checkpoint protein MAD2A; HsMAD2; Mitotic arrest deficient 2-like protein 1; MAD2-like protein 1
Human Swiss Prot No.:	Q13257
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAD2 Polyclonal Antibody detects endogenous levels of MAD2 protein.
Storage and Stability:	-20°C/1 year