

Antibodies Product Datasheet

Product Name:	GP-39 Antibody
Catalog Number:	C21661
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GP-39
Concentration:	1 mg/ml
Gene Name:	CHI3L1
Protein Name:	GP-39
Other Name:	Chitinase-3-like protein 1 (39 kDa synovial protein) (Cartilage glycoprotein 39) (CGP-39) (GP-39) (hCGP-39) (YKL-40)
Human Swiss Prot No.:	P36222
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GP-39
Storage and Stability:	-20°C/1 year