

Antibodies Product Datasheet

Product Name:	CD269 Antibody
Catalog Number:	C21673
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CD269.
Concentration:	1 mg/ml
Gene Name:	TNFRSF17
Protein Name:	Tumor necrosis factor receptor superfamily member 17
Other Name:	TNFRSF17; BCM; BCMA; Tumor necrosis factor receptor superfamily member 17; B-cell maturation protein; CD269
Human Swiss Prot No.:	Q02223
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD269 Polyclonal Antibody detects endogenous levels of CD269 protein.
Storage and Stability:	-20°C/1 year