

Antibodies Product Datasheet

Product Name:	KLK15 Antibody
Catalog Number:	C21951
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLK15
Concentration:	1 mg/ml
Gene Name:	KLK15
Protein Name:	KLK15
Other Name:	Kallikrein-15 (EC 3.4.21.-) (ACO protease)
Human Swiss Prot No.:	Q9H2R5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KLK15
Storage and Stability:	-20°C/1 year