

Antibodies Product Datasheet

Product Name:	Hamartin Antibody
Catalog Number:	C22116
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Hamartin . at AA range: 360-440
Concentration:	1 mg/ml
Gene Name:	TSC1 KIAA0243 TSC
Protein Name:	Hamartin
Other Name:	tuberous sclerosis 1
Human Swiss Prot No.:	Q92574
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Hamartin Polyclonal Antibody detects endogenous levels of Hamartin
Storage and Stability:	-20°C/1 year