

Antibodies Product Datasheet

Product Name:	Latexin Antibody
Catalog Number:	C22138
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Latexin . at AA range: 90-170
Concentration:	1 mg/ml
Gene Name:	LXN
Protein Name:	Latexin
Other Name:	LXN; Latexin; Endogenous carboxypeptidase inhibitor; ECI; Protein MUM; Tissue carboxypeptidase inhibitor; TCI
Human Swiss Prot No.:	Q9BS40
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Latexin Polyclonal Antibody detects endogenous levels of Latexin protein.
Storage and Stability:	-20°C/1 year