

Antibodies Product Datasheet

Product Name:	DR3 Antibody
Catalog Number:	C22143
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from DR3 . at AA range: 230-310
Concentration:	1 mg/ml
Gene Name:	TNFRSF25
Protein Name:	Tumor necrosis factor receptor superfamily member 25
Other Name:	TNFRSF25; APO3; DDR3; DR3; TNFRSF12; WSL; WSL1; Tumor necrosis factor receptor superfamily member 25; Apo-3; Apoptosis-inducing receptor AIR; Apoptosis-mediating receptor DR3; Apoptosis-mediating receptor TRAMP; Death receptor 3; Lymphocyte-associated rece
Human Swiss Prot No.:	Q93038
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DR3 Polyclonal Antibody detects endogenous levels of DR3 protein.
Storage and Stability:	-20°C/1 year