

Antibodies Product Datasheet

Product Name:	CTGF Antibody
Catalog Number:	C22169
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTGF
Concentration:	1 mg/ml
Gene Name:	CTGF CCN2 HCS24 IGFBP8
Protein Name:	CTGF
Other Name:	Connective tissue growth factor (CCN family member 2) (Hypertrophic chondrocyte-specific protein 24) (Insulin-like growth factor-binding protein 8) (IBP-8) (IGF-binding protein 8) (IGFBP-8)
Human Swiss Prot No.:	P29279
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CTGF
Storage and Stability:	-20°C/1 year