

Antibodies Product Datasheet

Product Name:	HSPBAP1 Antibody
Catalog Number:	C30002
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human HSPBAP1. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	HSPBAP1
Protein Name:	HSPB1-associated protein 1
Other Name:	HSPBAP1; PASS1; HSPB1-associated protein 1; 27 kDa heat shock protein-associated protein 1; Protein associated with small stress protein 1
Human Swiss Prot No.:	Q96EW2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSPBAP1 Polyclonal Antibody detects endogenous levels of HSPBAP1 protein.
Storage and Stability:	-20°C/1 year