

Antibodies Product Datasheet

Product Name:	LXR α Antibody
Catalog Number:	C30010
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human NR1H3. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	NR1H3
Protein Name:	Oxysterols receptor LXR-alpha
Other Name:	NR1H3; LXRA; Oxysterols receptor LXR-alpha; Liver X receptor alpha; Nuclear receptor subfamily 1 group H member 3
Human Swiss Prot No.:	Q13133
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LXR α Polyclonal Antibody detects endogenous levels of LXR α protein.
Storage and Stability:	-20°C/1 year