

Antibodies Product Datasheet

Product Name:	AChR α 9 Antibody
Catalog Number:	C30011
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CHRNA9. AA range:50-100
Concentration:	1 mg/ml
Gene Name:	CHRNA9
Protein Name:	cholinergic receptor, nicotinic, alpha 9
Other Name:	Neuronal acetylcholine receptor subunit alpha-9 (Nicotinic acetylcholine receptor subunit alpha-9) (NACHR alpha-9)
Human Swiss Prot No.:	Q9UGM1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AChR α 9 Polyclonal Antibody detects endogenous levels of cholinergic receptor, nicotinic, alpha 9
Storage and Stability:	-20°C/1 year