

Antibodies Product Datasheet

Product Name:	Cyclophilin D Antibody
Catalog Number:	C30012
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human PPID. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	PPID
Protein Name:	Peptidyl-prolyl cis-trans isomerase D
Other Name:	PPID; CYP40; CYPD; Peptidyl-prolyl cis-trans isomerase D; PPIase D; 40 kDa peptidyl-prolyl cis-trans isomerase; Cyclophilin-40; CYP-40; Cyclophilin-related protein; Rotamase D
Human Swiss Prot No.:	Q08752
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclophilin D Polyclonal Antibody detects endogenous levels of Cyclophilin D protein.
Storage and Stability:	-20°C/1 year