

Antibodies Product Datasheet

Product Name:	ET-1 Antibody
Catalog Number:	C30019
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human EDN1. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	EDN1
Protein Name:	Endothelin-1
Other Name:	EDN1; Endothelin-1; Preproendothelin-1; PPET1
Human Swiss Prot No.:	P05305
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ET-1 Polyclonal Antibody detects endogenous levels of ET-1 protein.
Storage and Stability:	-20°C/1 year