

Antibodies Product Datasheet

Product Name:	Eotaxin-3 Antibody
Catalog Number:	C30025
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CCL26. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	CCL26
Protein Name:	C-C motif chemokine 26
Other Name:	CCL26; SCYA26; C-C motif chemokine 26; CC chemokine IMAC; Eotaxin-3; Macrophage inflammatory protein 4-alpha; MIP-4-alpha; Small-inducible cytokine A26; Thymic stroma chemokine-1; TSC-1
Human Swiss Prot No.:	Q9Y258
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Eotaxin-3 Polyclonal Antibody detects endogenous levels of Eotaxin-3 protein.
Storage and Stability:	-20°C/1 year