

Antibodies Product Datasheet

Product Name:	IL-8 Antibody
Catalog Number:	C30028
Reactivity:	Human
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human IL8. AA range:50-99
Concentration:	1 mg/ml
Gene Name:	IL8
Protein Name:	Interleukin-8
Other Name:	IL8; CXCL8; Interleukin-8; IL-8; C-X-C motif chemokine 8; Emotakin; Granulocyte chemotactic protein 1; GCP-1; Monocyte-derived neutrophil chemotactic factor; MDNCF; Monocyte-derived neutrophil-activating peptide; MONAP; Neutrophil-activating protein 1; NA
Human Swiss Prot No.:	P10145
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-8 Polyclonal Antibody detects endogenous levels of IL-8 protein.
Storage and Stability:	-20°C/1 year