

Antibodies Product Datasheet

Product Name:	Bag-3 Antibody
Catalog Number:	C30035
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human BAG3. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	BAG3
Protein Name:	BAG family molecular chaperone regulator 3
Other Name:	BAG3; BIS; BAG family molecular chaperone regulator 3; BAG-3; Bcl-2-associated athanogene 3; Bcl-2-binding protein Bis; Docking protein CAIR-1
Human Swiss Prot No.:	O95817
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Bag-3 Polyclonal Antibody detects endogenous levels of Bag-3 protein.
Storage and Stability:	-20°C/1 year