

Antibodies Product Datasheet

Product Name:	FGF-12 Antibody
Catalog Number:	C30048
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FGF12. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	FGF12
Protein Name:	Fibroblast growth factor 12
Other Name:	FGF12; FGF12B; FHF1; Fibroblast growth factor 12; FGF-12; Fibroblast growth factor homologous factor 1; FHF-1; Myocyte-activating factor
Human Swiss Prot No.:	P61328
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-12 Polyclonal Antibody detects endogenous levels of FGF-12 protein.
Storage and Stability:	-20°C/1 year