

Antibodies Product Datasheet

Product Name:	Autotaxin Antibody
Catalog Number:	C30050
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human ENPP2. AA range:401-450
Concentration:	1 mg/ml
Gene Name:	ENPP2
Protein Name:	Ectonucleotide pyrophosphatase/phosphodiesterase family member 2
Other Name:	ENPP2; ATX; PDNP2; Ectonucleotide pyrophosphatase/phosphodiesterase family member 2; E-NPP 2; Autotaxin; Extracellular lysophospholipase D; LysoPLD
Human Swiss Prot No.:	Q13822
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Autotaxin Polyclonal Antibody detects endogenous levels of Autotaxin protein.
Storage and Stability:	-20°C/1 year