

Antibodies Product Datasheet

Product Name:	FGF-11 Antibody
Catalog Number:	C30053
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FGF11. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	FGF11
Protein Name:	Fibroblast growth factor 11
Other Name:	FGF11; FHF3; Fibroblast growth factor 11; FGF-11; Fibroblast growth factor homologous factor 3; FHF-3
Human Swiss Prot No.:	Q92914
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-11 Polyclonal Antibody detects endogenous levels of FGF-11 protein.
Storage and Stability:	-20°C/1 year