

## Antibodies Product Datasheet

|                               |                                                                                                                                             |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Annexin I Antibody                                                                                                                          |
| <b>Catalog Number:</b>        | C30064                                                                                                                                      |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                           |
| <b>Applications:</b>          | IF/ICC, WB, IHC-p, ELISA                                                                                                                    |
| <b>Dilutions:</b>             | IF: 1:50-200 WB 1:500-2000,IHC-p 1:500-200, ELISA 1:10000-20000                                                                             |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                            |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                           |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                       |
| <b>Immunogen:</b>             | Synthetic peptide from human protein at AA range: 130-180                                                                                   |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                     |
| <b>Gene Name:</b>             | ANXA1 ANX1 LPC1                                                                                                                             |
| <b>Protein Name:</b>          | Annexin A1 (Annexin I) (Annexin-1) (Calpactin II) (Calpactin-2) (Chromobindin-9) (Lipocortin I) (Phospholipase A2 inhibitory protein) (p35) |
| <b>Other Name:</b>            | Annexin A1 (Annexin I) (Annexin-1) (Calpactin II) (Calpactin-2) (Chromobindin-9) (Lipocortin I) (Phospholipase A2 inhibitory protein) (p35) |
| <b>Human Swiss Prot No.:</b>  | P04083                                                                                                                                      |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                     |
| <b>Specificity:</b>           | The antibody detects endogenous Annexin I                                                                                                   |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                |