

Antibodies Product Datasheet

Product Name:	GPR87/95 Antibody
Catalog Number:	C30073
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human GPR87. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	GPR87
Protein Name:	G-protein coupled receptor 87
Other Name:	GPR87; GPR95; FKSG78; G-protein coupled receptor 87; G-protein coupled receptor 95
Human Swiss Prot No.:	Q9BY21
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR87/95 Polyclonal Antibody detects endogenous levels of GPR87/95 protein.
Storage and Stability:	-20°C/1 year