

Antibodies Product Datasheet

Product Name:	FGF-4 Antibody
Catalog Number:	C30077
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FGF4. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	FGF4
Protein Name:	Fibroblast growth factor 4
Other Name:	FGF4; HST; HSTF1; KS3; Fibroblast growth factor 4; FGF-4; Heparin secretory-transforming protein 1; HST; HST-1; HSTF-1; Heparin-binding growth factor 4; HBGF-4; Transforming protein KS3
Human Swiss Prot No.:	P08620
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-4 Polyclonal Antibody detects endogenous levels of FGF-4 protein.
Storage and Stability:	-20°C/1 year